

**VENOSAN®**

**4000**

Compression Stockings

ADDITIONAL  
COLOURS  
FOR 2020

CLIMATE EFFECT



**VENOSAN®**  
Compression Stockings

[www.venosan.com](http://www.venosan.com)

VENOSAN® is a registered  
trademark of SWISSLASTIC AG  
ST. GALLEN, Switzerland

an  Company



SWISSLASTIC AG  
ST. GALLEN  
Rorschacher Strasse 304  
CH-9016 St. Gallen  
Switzerland



Swiss Made

**VENOSAN®**

# 4000

## Compression Stockings


From March 2020, VENOSAN® is expanding its VENOSAN® 4000 range with two additional colours. You can choose the stocking that best suits your individual style from a palette of six colours. The proven microfibre knit in the VENOSAN® 4000 compression stocking makes for lightness and well-being and from now on also in Marine and Anthracite.

VENOSAN® 4000 Compression Stockings	COMPRESSION CLASS	COLOUR	AD	AGH
	Ccl.II	Anthracite	▲	▲
		Marine	▲	▲

▲ only OTS (off the shelf)

### Colours

Due to technical printing reasons the colours may vary.



71 % TACTEL®  
Polyamide  
29 % Elastane




TACTEL® fibre is a registered trademark of INVISTA and is used under license to NILIT® for nylon fibre products.

# 4002

## Compression Stockings Ccl.2

VENOSAN® 4002  
Ccl.II  
23 – 32 mm Hg  
(3,1 – 4,3 kPa)



ARTICLE	COLOUR	LENGTH	TOE	SELF SUPPORTING TOP	S	M	L	XL	
AD	Anthracite	long	open		SA 42001	SA 42002	SA 42003	SA 42004	
		short			SA 42021	SA 42022	SA 42023	SA 42024	
		long	closed		SA 42011	SA 42012	SA 42013	SA 42014	
		short			SA 42031	SA 42032	SA 42033	SA 42034	
	Marine	long	open		SI 42001	SI 42002	SI 42003	SI 42004	
		short			SI 42021	SI 42022	SI 42023	SI 42024	
		long	closed		SI 42011	SI 42012	SI 42013	SI 42014	
		short			SI 42031	SI 42032	SI 42033	SI 42034	
AGH	Anthracite	long	open	plain	SA 42611	SA 42612	SA 42613	SA 42614	
		short			SA 42641	SA 42642	SA 42643	SA 42644	
		long	closed	plain	SA 42621	SA 42622	SA 42623	SA 42624	
		short			SA 42651	SA 42652	SA 42653	SA 42654	
		long		lace	SA 42601	SA 42602	SA 42603	SA 42604	
		short			SA 42631	SA 42632	SA 42633	SA 42634	
		Marine	long	open	plain	SI 42611	SI 42612	SI 42613	SI 42614
			short			SI 42641	SI 42642	SI 42643	SI 42644
	long		closed	plain	SI 42621	SI 42622	SI 42623	SI 42624	
	short				SI 42651	SI 42652	SI 42653	SI 42654	
	long			lace	SI 42601	SI 42602	SI 42603	SI 42604	
	short				SI 42631	SI 42632	SI 42633	SI 42634	